## RAPID SARS-COV-2 ANTIGEN TEST CARD

INSTRUCTION GUIDE FOR ANTERIOR NASAL SWAB SPECIMENS

For self-testing  
REF 1N40C5-8  
REF 1N40C5-9  
REF 1N40C5-10

CE  
0123  
IVD

Please follow the instruction leaflet carefully.

**INTENDED USE**  
Rapid SARS-CoV-2 Antigen Test Card is an immunochromatography based one step in vitro test. It is designed for the rapid qualitative determination of SARS-CoV-2 antigen in anterior nasal swabs from individuals with or without symptoms or other epidemiological reasons to suspect COVID-19. Rapid SARS-CoV-2 Antigen Test Card detects the SARS-CoV-2 nucleocapsid protein (N protein). Rapid SARS-CoV-2 Antigen Test Card shall not be used as sole basis to diagnose or exclude SARS-CoV-2 infection. Children under 14 years of age should be assisted by an adult.

**SUMMARY**  
The novel coronaviruses belong to the *β* genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronaviruses are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

**MATERIALS PROVIDED**

Components	For 1 Test/Box	For 5 Tests/Box	For 20 Tests/Box
Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)	1	5	20
Sterile swab	1	5	20
Extraction buffer tube	1	5	20
Instructions for use (this leaflet)	1	1	1
Tube stand	1 (packaging)	1	1

**PERFORMANCES (SENSITIVITY AND SPECIFICITY)**  
Rapid SARS-CoV-2 Antigen Test Card was compared to the confirmed clinical diagnosis. The study involved 1063 nasal samples. The test results are summarized below:

Evaluated Reagent Results	RT-PCR Results	Total
Positive (+)	425	1
Negative (-)	10	627
Total	435	628

Sensitivity (PPA) =  $425/435 \times 100\% = 97.70\% (95\% \text{ CI: } 96.29\% - 99.11\%)$

Specificity (NPV) =  $627/628 \times 100\% = 99.84\% (95\% \text{ CI: } 99.53\% - 99.99\%)$

Accuracy (OPA) =  $1052/1063 \times 100\% = 98.97\% (95\% \text{ CI: } 98.36\% - 99.57\%)$

A negative result does not exclude COVID-19.

- 99.84% of non-professionals carried out the test without requiring assistance.

- 99.92% of the different types of results were interpreted correctly.

**INTERFERENCES**

None of the following substances at the tested concentration showed any interference with the test.

Whole Blood: 1% Alkalot: 10% Mucin: 2% Phenylethanol: 15% Tobramycin: 0.0004% Octyleneglycin: 15% Benzalkonium: 0.15% Cetacrylate: 5% Mupirocin: 0.25% Zicam Nasal Spray: 5% Osselamivir Phosphate: 0.5% Sodium chloride: 5% Human Anti-mouse Antibody (HAMA): 60 ng/mL Biotin: 1200 ng/mL

**CROSS REACTIVITY**

Test Device was evaluated by testing viruses and other microorganisms. The final test concentrations of viruses and other microorganisms are documented in the Cross-Reactivity Study. The following viruses and other microorganisms except the Human coronavirus have no effect on the test results of the Test Device: Human coronavirus 229E, Human coronavirus OC43, Human coronavirus NL63, Human coronavirus HKU1, MERS coronavirus, Parainfluenza Virus 1-4, Enterovirus EV71, Respiratory syncytial virus, Rhinovirus, Influenza A virus (H1N1 and H3N2), Influenza B virus (Yamagata and Victoria), Adeno virus 71, Human Metapneumovirus (hMPV), Staphylococcus epidemidis, Chlamydia pneumoniae, Mycoplasma pneumoniae, Streptococcus pyogenes, Bordetella pertussis, Mycobacterium tuberculosis, Legionella pneumonia, Mycoplasma pneumoniae, Haemophilus influenzae, Candida albicans, Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, Pneumocystis jirovecii (PJP) and Pneumocystis jirovecii (PJP).

**IMPORTANT INFORMATION BEFORE THE EXECUTION**

1. Read this instruction guide carefully.

2. Do not use the product beyond the expiration date.

3. Do not use the product if the seal is damaged or broken.

4. Do not use the test device at 4 to 10 °C until the original sealed pouch. Do Not Freeze.

5. The product should be used at room temperature (15°C to 30°C). If the product has been stored in a cool area (less than 15°C), leave it at normal room temperature for 30 minutes before using.

6. Handle all specimens as potentially infectious.

7. Do not use the test device for specimen collection, storage, and transport may yield inaccurate test results.

8. Correct specimen collection is the most important step in the procedure. Make sure to collect enough specimen material (nasal secretion) with the swab, especially for anterior nasal sampling.

9. Blow the nose several times before collecting specimen.

10. Take care when handling the test strip to prevent damage to the specimen well (S).

11. Apply the drops of test solution only to the specimen well (S).

12. Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.

13. When used as intended, there should not be any contact with the extraction buffer. In case of contact with skin, eyes, mouth or other parts, rinse with clean water. If an irritation persists, consult a medical professional.

14. Children under 14 years of age should be assisted by an adult.

**LIMITATIONS**

1. The test is to be used exclusively for the qualitative detection of SARS-CoV-2 viral antigen in anterior nasal swab specimens. The exact concentration of SARS-CoV-2 viral antigen cannot be determined as part of this test.

2. Proper specimen collection is critical. Failure to follow the procedure may result in inaccurate test results. Improper collection, storage, or even freezing and melting of the specimen will affect the test results.

3. If the viral load of the specimen is below the detection limit of the test, the test may produce a negative result.

4. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be made by the physician after evaluation of all clinical and laboratory results.

5. A negative result does not exclude viral infection except for SARS-CoV-2 and should be confirmed by molecular diagnostic methods if COVID-19 is suspected.

6. A positive result does not exclude co-infection with other pathogens.

7. The SARS-CoV-2 rapid antigen test can detect both viable and non-viable SARS-CoV-2 material. The performance of the SARS-CoV-2 rapid test is dependent on viral load and may correlate with other diagnostic methods performed on the same specimen.

8. Use the test device within 10 minutes of specimen collection.

9. Sensitivity for nasal or oropharyngeal swabs may be lower than nasopharyngeal swabs. It is recommended to use the nasopharyngeal swab specimens by healthcare professionals.

10. Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope.

11. The amount of antigen in a sample may decrease as the duration of illness increases. Specimens collected after day 5-7 of illness are more likely to be tested negative compared to a RT-PCR assay.

12. The kit was validated with the assayed swabs. Use of alternative swabs may result in false negative results.

13. The ability of Rapid SARS-CoV-2 Antigen Test Card has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.

14. Positive test results do not rule out co-infections with other pathogens. Positive results may occur in cases of infection with SARS-CoV.

**PREPARATION**

● Clear, clean and dry a flat surface.

● Check the test kit contents. Make sure that nothing is damaged or broken.

● Timer at hand.

● Wash your hands.

**DISPENSAL**


The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

**PROCEDURE**

This test is suitable for people of all ages. The recommended operators are aging from 14 to 90. Children under 14 years of age should be tested by an adult. Do not continue the test if the child feels any pain.


1 Remove the sealing from the extraction buffer tube.

Caution: Open it away from your face and be careful not to spill any of the liquid.

2 Find the swab in the sealed wrapper in front of you. Identify the soft, fabric tip of the swab.

Caution: Never touch the soft, fabric tip of the swab with your hands.


3 Carefully insert swab into one nostril. The swab tip should be inserted no less than 2.5 cm (1 inch) from the edge of the nostril. Roll swab 3-4 times along the mucosa inside the nostril. Leave swab in the nostril for several seconds. Using the same swab, repeat this process for the other nostril. Withdraw swab from the nasal cavity.

Caution: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.


4 Place the swab with specimen into the extraction tube. Roll the swab three to five (3-5) times. Leave swab in the extraction buffer for 1 minute.

5 Pinch extraction tube with fingers and remove the solution from swab as much as possible.

6 Install the nozzle cap onto the extraction tube lightly.


7 Bring the kit components to room temperature before testing. Open the pouch and remove the card. Place the card on a flat and level surface.

Caution: Once opened, the test card must be used immediately.

8 Invert the extraction tube and add 3 drops (about 75 µl) of test specimen into the specimen well (S) by gently squeezing the extraction tube.


Caution: The formation of air bubbles in the specimen well (S) must be avoided. If there are air bubbles, add another drop of test specimen.

9 Read the results at 15-20 minutes.

Caution: Results after 20 minutes may not be accurate.

The used device may be disposed of with normal household waste in accordance with the applicable local regulations.


**INTERPRETATION OF RESULTS**

**Positive:**


If two colored bands appear with one colored band in the Control Zone (C) and another in the Test Zone (T) within 15-20 minutes, the test result is positive.

Caution: No matter how faint the colored band is in the Test Zone (T), the result should be considered as positive.


**Negative:**

If no color line appears in the control area (C) within 15-20 minutes, the test is invalid. Repeat the test with a new test card.

**Invalid:**

If no color line appears in the control area (C) within 15-20 minutes, the test is invalid. Repeat the test with a new test card.

**QUALITY CONTROL**

The control line is an integrated reagent and is used to control the procedure. The control line appears when the test has been performed correctly and the reagents are reactive.

**FREQUENTLY ASKED QUESTIONS (FAQ)**

1. How does the detection work?

The N protein of SARS-CoV-2 virus reacts with the stripe-like coating of the test line and, if present, results in a color change, i.e. a red line appears.

2. When should I test myself?

You can test yourself whether you have symptoms or not. Studies show that earlier testing within the first 4 days of illness typically means a higher viral load, which is easier to detect. Since the test result is a snapshot valid for that point in time, testing should be repeated as recommended by local authorities.

3. How often should I test myself? What should I do?

Be sure to follow the specific multiple times before collecting the specimen.

4. The test strip is clearly discolored or smudged? What is the reason for this?

Please note that the test card should not be used if the test strip shows the liquid absorption of the test strip is naturally limited. If the control line does not appear, the test strip is discolored or dislodged, making it unreadable, please repeat the test according to the instructions.

5. I have taken the test, but I don't see a control line (C). What should I do?

Your test result is invalid. Observe the answer to question 4 and repeat the test according to the instructions for use.

6. I am unsure about reading the result. What should I do?

For the test to be positive, the two lines must be clearly visible with the full width of the cassette. If you are still unsure about the results, contact the nearest medical facility according to the recommendations of your local authorities.

7. My result is positive. What should I do?

If your result is positive and the test strip clearly indicates the control line as well as the test line, you should contact the nearest medical facility as recommended by your local authorities. Your test result may be double-checked and the authority or facility will explain the appropriate next steps.

8. My result is negative. What should I do?

If the test kit only clearly shows the control line, this may mean that you are negative or that the viral load

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